



### Intelli-Temp Calibrator

- Calibrate the SMART System<sup>5</sup> temperature controller to maintain optimum performance
- Easy-to-use
- NIST traceable
- ISO procedures
- Fine-tune your temperature feedback system in minutes!

The Intelli-Temp Calibrator heats to precisely 130 °C and maintains that temperature while the integrated software calibrates the temperature sensor in the analyzer. It takes only minutes to keep your SMART System in the best possible operating condition.



### Bar Code Reader

- High resolution scanner for UPC Codes
- Scan sample identification numbers quickly and accurately
- Ideal for laboratories with a heavy workflow

Scan UPC Codes for samples quickly and easily with the handheld Bar Code Reader. Integrated software stores sample numbers and test results.



*Innovators in  
Microwave Technology*

(800) 726-3331 • (704) 821-7015 • Fax (704) 821-5185  
www.cem.com • e-mail: info@cem.com

